

User Setup		Vault Permissions		Permissions on Assigned Study							Additional Information
Security Profile	Application Role on Study	Manage Vault Users	Access Vault Reports	Manage Study Users	Manage Study	Manage Study Country	Manage Site	Access Studio	Access Study Home	Lock and Unlock Sites and Studies	Notes
Vault Owner	Not applicable	Yes	Edit	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Users with the Vault Owner security profile have all permissions on all studies.
System Administrator	Not applicable	Yes	Edit	Yes	Yes	Yes	Yes	No	No	No	None
Clinical Administrator	If not assigned to study	No	Edit	Yes	Yes	Yes	Yes	No	No	No	This is the typical configuration for a central administration function.
	Monitor	No	Edit	Yes	Yes	Yes	Yes	No	Yes	No	None
	Study Manager	No	Edit	Yes	Yes	Yes	Yes	No	Yes	No	None
	Study Viewer	No	Edit	Yes	Yes	Yes	Yes	No	Yes	No	We do not recommend this combination because a clinical administrator will always be able to manage study country and site, which is more permission than a study viewer is expected to have.
	Data Manager	No	Edit	Yes	Yes	Yes	Yes	Yes	Yes	No	A user with this permission has write access in Studio for Exports. They have read-only access to Collections and all other Study Settings.
	Lead Data Manager	No	Edit	Yes	Yes	Yes	Yes	Yes	Yes	Yes	A user with this permission has write access to Study Settings - Deployments, and Exports. They have read-only access to Collections and all other Study Settings.
Study Team User	Study Builder	No	Edit	Yes	Yes	Yes	Yes	Yes	No	Yes	None
	If not assigned to study	No	View/Run	No	No	No	No	No	No	No	A user with the Read-only User security profile has no permissions on this study until assigned to an application role.
	Data Manager	No	View/Run	Yes	No	No	No	Yes	Yes	No	A user with this permission has write access in Studio for Exports. They have read-only access to Collections and all other Study Settings.
	Lead Data Manager	No	View/Run	Yes	No	No	No	Yes	Yes	Yes	A user with this permission has write access to Study Settings - Deployments, and Exports. They have read-only access to Collections and all other Study Settings.
	Monitor	No	View/Run	Yes	Yes	Yes	Yes	No	Yes	No	None
	Study Builder	No	View/Run	Yes	Yes	No	No	Yes	No	Yes	None
	Study Manager	No	View/Run	Yes	Yes	Yes	Yes	No	Yes	No	None
Study Viewer	No	View/Run	No	No	No	No	No	Yes	No	None	
Read-only User	If not assigned to study	No	View/Run	No	No	No	No	No	No	No	A user with the Read-only User security profile has no permissions on this study until assigned to an application role.
	Data Manager	No	View/Run	No	No	No	No	No	Yes	No	None
	Monitor Study Manager Study Viewer Lead Data Manager	No	View/Run	No	No	No	No	No	Yes	No	None
	Study Builder	No	View/Run	No	No	No	No	No	No	No	We do not recommend this combination because a user